



# Chemistry Analyzer





## On-Site Diagnostic Testing, with 10 Panels to Choose From


The VetScan VS2 is a state-of-the-art chemistry, electrolyte, immunoassay and blood gas analyzer that delivers uncompromising accuracy from just two drops of whole blood, serum or plasma. The VetScan VS2 is amazingly simple and intuitive, with precision reference laboratory quality results in 12 minutes—making it ideal for veterinary clinics, mobile practitioners, research laboratories, and pharmaceutical and biotech companies.

## Features

- Ease of use
- Expansive and unique menu
- Accuracy comparable to reference laboratories
- Intelligent Quality Control (iQC)
- Little to no maintenance
- Detects sample interference - hemolysis, lipemia, and icterus
- Small footprint and portable (only 11.6 lbs.)
- Cost effective
- Speed: 12 minutes to results
- Stores 5,000 patient records
- Connectivity options

## Benefits

- Reduce training time and increase employee compliance to use
- Save time, counter space, and costs with ability to test chemistry and electrolyte analytes on a single panel and have access to specialty tests including T<sub>4</sub>, Bile Acids, TCO<sub>2</sub>, and Canine Heartworm on one instrument
- Diagnose and treat patients confidently
- No more worrying about running controls and calibrations as iQC monitors the analyzer, rotor, and sample before and during every analysis to ensure reliable results
- Save time and costs plus increase employee compliance to use
- Diagnose and treat patients confidently as analyte results are suppressed if affected by > 10% interference and receive useful information about sample integrity and handling
- Use this analyzer anywhere and anytime in various environments
- Save money with low-price reagents, no control and calibration costs, and no service contracts
- Obtain results, diagnose, discuss with clients and begin treatment within one visit
- Provides peace of mind that results are secure and allows recall of patient history
- Connect with Practice Management Software and external printers using the four USB ports and LAN port



# Comprehensive Test Portfolio – 10 Different Profiles to Choose From

Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	CL-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHOS	T4	TBIL	tCO <sub>2</sub>	TP	UA
Canine Wellness Profile** #500-0044	•	•	•	/	/	/	•	•	/	•	/	/	•	/	•	•	/	/	/	•	/	•	/	•	/
Comprehensive Diagnostic Profile #500-0038	•	•	•	•	/	/	•	•	/	/	/	/	•	/	•	•	•	/	•	•	/	•	/	•	/
Prep Profile II #500-0026	/	•	•	/	/	/	•	/	/	/	/	/	•	/	/	•	/	/	/	/	/	/	/	•	/
Critical Care Plus #500-0042	/	/	•	/	/	/	•	/	/	/	/	•	•	/	/	•	•	/	•	/	/	/	•	/	/
T4/Cholesterol Profile #500-0037	/	/	/	/	/	/	/	/	•	/	/	/	/	/	/	/	/	/	/	/	•	/	/	/	/
Kidney Profile Plus** #500-0045	•	/	/	/	/	/	•	•	/	/	/	•	•	/	/	•	•	/	•	•	/	/	•	/	/
Mammalian Liver Profile #500-0040	•	•	•	/	/	•	•	/	•	/	/	/	/	•	/	/	/	/	/	/	/	/	•	/	/
Equine Profile Plus #500-0043	•	/	/	/	•	/	•	•	/	/	•	/	•	•	•	•	•	/	•	/	/	•	•	•	/
Avian/Reptilian Profile Plus #500-0041	•	/	/	/	•	•	/	•	/	/	•	/	/	/	•	•	•	/	•	•	/	/	/	•	•
Large Animal Profile #500-0023	•	•	/	/	•	/	•	•	/	/	•	/	/	•	•	/	/	•	/	•	/	/	/	•	/

\*Calculated value

\*\*For VetScan VS2 Only

## Intelligent Quality Control (iQC)

With the VetScan VS2 Chemistry Analyzer, quality control is built-in—thanks to a sophisticated iQC system that combines with proprietary algorithms to assure quality and dependability with every result. No more worrying about running controls and calibrations.



## Simple Operation


Practices everywhere have discovered how the VetScan Chemistry Analyzer has improved practice efficiencies and saved them time and money. Simplicity, speed, through-put with results on-site, in minutes.



1 Insert sample.



2 Insert rotor.



3 Read results.

Requires only two drops of whole blood.



## VetScan VS2

### Easy as 1-2-3

For more information on the VetScan VS2, contact your Abaxis Representative at 800.822.2947, email us at [vetscan@abaxis.com](mailto:vetscan@abaxis.com) or click [www.abaxis.com](http://www.abaxis.com).

Abaxis is your total on-site laboratory partner, offering the comprehensive and innovative VetScan line of diagnostic instruments and individual tests. The VetScan name has become synonymous with convenience, ease of use and cost-effectiveness, delivering on our promise of speed, accuracy and reliability. It's not just better technology. It's a better way.

For questions:

 800.822.2947  [www.abaxis.com](http://www.abaxis.com)  [vetscan@abaxis.com](mailto:vetscan@abaxis.com)



#### Abaxis, Inc.

3240 Whipple Road  
Union City, CA 94587  
Tel 800.822.2947  
Fax 510.441.6150  
Email [vetscan@abaxis.com](mailto:vetscan@abaxis.com)

[www.abaxis.com](http://www.abaxis.com)

#### ABAXIS Europe GmbH

Bunsenstr. 9-11  
64347 Griesheim Germany  
Tel +49 6155 780 210  
Fax +49 6155 780 2111  
Email [abaxis@abaxis.de](mailto:abaxis@abaxis.de)

GLOBAL DIAGNOSTICS

